



LX9 PNEUMATIC COMPRESSION PUMPS

- Reduces limb volume and size
- Softens tissue and improves skin quality
- Cost effective, portable, light and comfortable to wear
- A great treatment option for home management





SEQUENTIAL PNEUMATIC COMPRESSION PUMPS

Intermittent pneumatic compression mimics the body's manual lymphatic drainage. By applying pressure to the surface of the limb it promotes the circulation of blood/fluid throughout the body. Each chamber inflates and deflates sequentially starting from the extremity and moving towards the centre of the body.

LX9 FEATURES

- Variable pressure 10-180mmHg
- Continuous mode - 80mmHg max pressure
- Therapy timer & lockable programs
- Easy fit garments - leg, arm or trunk
- Various sizes
- Four treatment modes
- Six cycle speeds
- Sanitary, washable boot-liners available

RENT-TRY-BUY

We offer a Rent Try Buy scheme.

If you decide to purchase, we will refund up to three months of paid rental off the purchase price.

Australia

☎ 1300 726 666

✉ info@medirent.com.au

🌐 www.medirent.com.au

New Zealand

☎ 0800 725 555

✉ info@medirent.com.au

🌐 www.medirent.co.nz